
PRAMA


Product Overview



FUTURE
DENTAL

Neck morphology

Prama, the one and only intramucosal implant, was designed with a converging neck to maximize the thickness of soft tissues. To answer all clinical situations, Prama implants are available in three different heights of the intramucosal neck:




Short Neck
neck h 1.80 mm

Regular Neck
neck h 2.80 mm

Long Neck
neck h 3.80 mm

Prama Slim

Prama Slim was specifically developed for those areas with limited horizontal dimension. This implant has a straight intramucosal neck 1.80 mm or 2.80 mm high.






ø 3.30 Short Neck
neck h 1.80 mm

ø 3.30 Regular Neck
neck h 2.80 mm

Note: The use of Prama Slim for single crowns is just intended for the upper lateral incisors and for the lower central and lateral incisors. In all the other positions it can only be used for multiple rehabilitations.

Neck morphology and prosthesis

The concept of intramucosal Implant is very simple from the prosthetic point of view: considering the soft tissues thickness, it will be chosen the height of the neck between 1.80 mm, 2.80 mm and 3.80 mm and it will be decided whether to use a prosthesis that lays on the implant platform or whether to partially embrace the converging neck with dedicated Prama IN components.
All three necks can receive the prosthesis that closes at the platform level, while the Prama IN prosthesis that embraces the neck for 0.50 mm can only be used with 2.80 mm and 3.80 mm necks.
The Prama IN components that embrace the neck for 1.50 mm in height can only be used on the 2.80 mm neck.

		Standard prosthesis laying on the implant platform	0.50 mm Prama IN prosthesis	1.50 mm Prama IN prosthesis
				
Prama Slim	Collo			
	1.80  ø 3.30	✓	✗	✗
	2.80  ø 3.30	✓	✗	✗
Prama, Prama RF and Prama RF SL	1.80  ø 3.80 ø 4.25	✓	✗	✗
	2.80  ø 3.80 ø 4.25 ø 5.00	✓	✓	✓
	3.80  ø 3.80 ø 4.25 ø 5.00	✓	✓	✗

Guide to implant codes

Prama with cylindrical morphology

Morphology and neck	Surface	Diameter	Length
LAS: 1.80 mm - Short Neck **LA**: 2.80 mm - Regular Neck **LAL**: 3.80 mm - Long Neck example: **LAS-ZT-380-100**	**ZT**: ZirTi surface and UTM neck	**330**: 3.30 mm - Prama Slim **380**: 3.80 mm **425**: 4.25 mm **500**: 5.00 mm	**060**: 6.00 mm **085**: 8.50 mm **100**: 10.00 mm **115:** 11.50 mm **130**: 13.00 mm **150**: 15.00 mm
LAS-	**ZT-**	**380-**	**100**

Prama with conical morphology

Morphology and neck	Surface	Diameter	Spira		Length
			Regular	Large	
LSS: 1.80 mm - Short Neck **LS**: 2.80 mm - Regular Neck **LSL**: 3.80 mm - Long Neck example: **LSL-ZT-380SL-100**	**ZT**: ZirTi surface and UTM neck	**380**: 3.80 mm **425**: 4.25 mm **500**: 5.00 mm	**-**	**SL**	**060**: 6.00 mm **085**: 8.50 mm **100**: 10.00 mm **115:** 11.50 mm **130**: 13.00 mm **150**: 15.00 mm
LSL-	**ZT-**	**380**		**SL-**	**100**

Prama Slim – Short and Regular Neck



Prama Slim
Short Neck

h	ø 3.30
6.00	
8.50	**LAS-ZT-330-085**
10.00	**LAS-ZT-330-100**
11.50	**LAS-ZT-330-115**
13.00	**LAS-ZT-330-130**
15.00	**LAS-ZT-330-150**

Surgical cover screw	
	L-VT-340



Prama Slim
Regular Neck

h	ø 3.30
6.00	**LA-ZT-330-060**
8.50	**LA-ZT-330-085**
10.00	**LA-ZT-330-100**
11.50	**LA-ZT-330-115**
13.00	**LA-ZT-330-130**
15.00	**LA-ZT-330-150**

Surgical cover screw	
	L-VT-340

Prama – Short Neck



Prama

h	ø 3.80	ø 4.25
6.00		
8.50	**LAS-ZT-380-085**	**LAS-ZT-425-085**
10.00	**LAS-ZT-380-100**	**LAS-ZT-425-100**
11.50	**LAS-ZT-380-115**	**LAS-ZT-425-115**
13.00	**LAS-ZT-380-130**	**LAS-ZT-425-130**
15.00	**LAS-ZT-380-150**	**LAS-ZT-425-150**

Surgical cover screw	
	L-VT-340



Prama RF

h	ø 3.80	ø 4.25
6.00		
8.50	**LSS-ZT-380-085**	**LSS-ZT-425-085**
10.00	**LSS-ZT-380-100**	**LSS-ZT-425-100**
11.50	**LSS-ZT-380-115**	**LSS-ZT-425-115**
13.00		
15.00		

Surgical cover screw	
	L-VT-340

Prama - Regular Neck



Prama

h	ø 3.80	ø 4.25	ø 5.00
6.00	LA-ZT-380-060	LA-ZT-425-060	LA-ZT-500-060
8.50	LA-ZT-380-085	LA-ZT-425-085	LA-ZT-500-085
10.00	LA-ZT-380-100	LA-ZT-425-100	LA-ZT-500-100
11.50	LA-ZT-380-115	LA-ZT-425-115	LA-ZT-500-115
13.00	LA-ZT-380-130	LA-ZT-425-130	LA-ZT-500-130
15.00	LA-ZT-380-150	LA-ZT-425-150	LA-ZT-500-150



Prama RF

h	ø 3.80	ø 4.25	ø 5.00
6.00		LS-ZT-425-060	LS-ZT-500-060
8.50	LS-ZT-380-085	LS-ZT-425-085	LS-ZT-500-085
10.00	LS-ZT-380-100	LS-ZT-425-100	LS-ZT-500-100
11.50	LS-ZT-380-115	LS-ZT-425-115	LS-ZT-500-115
13.00	LS-ZT-380-130	LS-ZT-425-130	LS-ZT-500-130
15.00	LS-ZT-380-150	LS-ZT-425-150	LS-ZT-500-150



Prama RF SL

h	ø 3.80	ø 4.25	ø 5.00
6.00			
8.50	LS-ZT-380SL-085	LS-ZT-425SL-085	LS-ZT-500SL-085
10.00	LS-ZT-380SL-100	LS-ZT-425SL-100	LS-ZT-500SL-100
11.50	LS-ZT-380SL-115	LS-ZT-425SL-115	LS-ZT-500SL-115
13.00	LS-ZT-380SL-130	LS-ZT-425SL-130	LS-ZT-500SL-130
15.00	LS-ZT-380SL-150	LS-ZT-425SL-150	LS-ZT-500SL-150

Surgical cover screw	
	L-VT-340

Prama - Long Neck



Prama

h	ø 3.80	ø 4.25	ø 5.00
6.00	LAL-ZT-380-060	LAL-ZT-425-060	LAL-ZT-500-060
8.50	LAL-ZT-380-085	LAL-ZT-425-085	LAL-ZT-500-085
10.00	LAL-ZT-380-100	LAL-ZT-425-100	LAL-ZT-500-100
11.50	LAL-ZT-380-115	LAL-ZT-425-115	LAL-ZT-500-115
13.00	LAL-ZT-380-130	LAL-ZT-425-130	
15.00			



Prama RF

h	ø 3.80	ø 4.25	ø 5.00
6.00		LSL-ZT-425-060	LSL-ZT-500-060
8.50	LSL-ZT-380-085	LSL-ZT-425-085	LSL-ZT-500-085
10.00	LSL-ZT-380-100	LSL-ZT-425-100	LSL-ZT-500-100
11.50	LSL-ZT-380-115	LSL-ZT-425-115	LSL-ZT-500-115
13.00	LSL-ZT-380-130	LSL-ZT-425-130	
15.00			



Prama RF SL

h	ø 3.80	ø 4.25	ø 5.00
6.00			
8.50	LSL-ZT-380SL-085	LSL-ZT-425SL-085	LSL-ZT-500SL-085
10.00	LSL-ZT-380SL-100	LSL-ZT-425SL-100	LSL-ZT-500SL-100
11.50	LSL-ZT-380SL-115	LSL-ZT-425SL-115	LSL-ZT-500SL-115
13.00			
15.00			

Surgical cover screw	
	L-VT-340

Surgical kits

Prama, Prama RF and Prama RF SL surgical kit



ZPRAMA-INT	complete surgical grommetless kit of the instruments necessary for Prama and Prama RF implants
L-TRAYL-INT	radel instrument grommetless tray

Drilling kit for Prama cylindrical short implants



ZSHORTY-INT	drilling kit Shorty for Prama short implants
SHORTY-TRAY-INT	radel instrument tray
GROMMET-CA-1	spare grommets for drills or contra-angle instrumentsts (5 units package)

Please note: the ZSHORTY-INT kit allows you to insert the Shorty implants of the Premium, Outlink² and Prama implant lines: please refer to the MC-IMP-PRAMA surgical manual for the correct sequence of use of the only instruments necessary for Prama diameters. ZSHORTY-INT kit does not allow the insertion of Prama RF short implants, for which the instruments contained in the L-INTEGRA-060 set at page 8 are necessary.


Surgical kits for guided insertion of Prama implants



ZGS-PRO-A-INT	grommetless surgical kit including all the instruments for guided insertion of Prama implants
GS-PROA-INTEGRA-060	set of drills for the preparation of Prama implants, height 6.00 mm
GSPROA-TRAY-INT	radel instrument grommetless tray

For further information on guided insertion of implants, please consult the surgical manual MC-UD-EPLANPROA-E.

Surgical kits for guided insertion of Prama RF and Prama RF SL implants



ZGS-PRO-S-INT	grommetless surgical kit including all the instruments for guided insertion of Prama RF and Prama RF SL implants
GS-PROS-INTEGRA-060	set of drills for the preparation of Prama RF implants, height 6.00 mm
GSPROS-TRAY-INT	radel instrument grommetless tray

For further information on guided insertion of implants, please consult the surgical manual MC-UD-EPLANPROS-E.

Instruments included in Prama surgical kit

Precision drill

FS-230

Pilot drill and related stops



FPT3-200-LXS



h	
8.50	**STOP4-200-085**
10.00	**STOP4-200-100**
11.50	**STOP4-200-115**
13.00	**STOP4-200-130**
15.00	**STOP4-200-150**

Intermediate drills and related stops



FFT3-250-LXS 2.50 cylindrical drill

SE-FK250 2.50 conical drill



h	cylindrical	conical
8.50	**STOP4-250-085**	**CSR-STOP-2028-085**
10.00	**STOP4-250-100**	**CSR-STOP-2028-100**
11.50	**STOP4-250-115**	**CSR-STOP-2028-115**
13.00	**STOP4-250-130**	**CSR-STOP-2028-130**
15.00	**STOP4-250-150**	**CSR-STOP-2028-150**

Prama RF Reply



h	ø 3.80	ø 4.25	ø 5.00
8.50	**SH-380-085-RP**	**SH-425-085-RP**	**SH-500-085-RP**
10.00	**SH-380-100-RP**	**SH-425-100-RP**	**SH-500-100-RP**
11.50	**SH-380-115-RP**	**SH-425-115-RP**	**SH-500-115-RP**
13.00	**SH-380-130-RP**	**SH-425-130-RP**	**SH-500-130-RP**
15.00	**SH-380-150-RP**	**SH-425-150-RP**	**SH-500-150-RP**

Prama bone taps with hexagonal connector



ø 3.30	ø 3.80	ø 4.25	ø 5.00
A-MS-330	**A-MS-380**	**A-MS-425**	**A-MS-500**

Prama RF bone taps with right angle shank



ø 3.80	ø 4.25	ø 5.00
SH-MS-380-CA	**SH-MS-425-CA**	**SH-MS-500-CA**

Screwdrivers



HSMXS-20-DG	digital, extra short
HSM-20-DG	digital, short
HSML-20-DG	digital, long



HSM-20-EX	with hexagonal connector
HSML-20-EX	with hexagonal connector, long



HSM-20-CA	with right angle shank

Dynamometric ratchet



CRI5-KIT

Pins for neck positioning



h	
8.50	**L-PP-085**
10.00	**L-PP-100**
11.50	**L-PP-115**
13.00	**L-PP-130**
15.00	**L-PP-150**

Parallelism pin with depth laser marking



PPTL-2-28

Prama final drills and related stops



ø 3.30	ø 3.80	ø 4.25	ø 5.00
FFT3-280-LXS	**FFT3-300-LXS**	**FFT3-340-LXS**	**FFT3-425-LXS**



h	ø 3.30	ø 3.80	ø 4.25	ø 5.00
8.50	**STOP4-280-085**	**STOP4-300-085**	**STOP4-340-085**	**STOP4-425-085**
10.00	**STOP4-280-100**	**STOP4-300-100**	**STOP4-340-100**	**STOP4-425-100**
11.50	**STOP4-280-115**	**STOP4-300-115**	**STOP4-340-115**	**STOP4-425-115**
13.00	**STOP4-280-130**	**STOP4-300-130**	**STOP4-340-130**	**STOP4-425-130**
15.00	**STOP4-280-150**	**STOP4-300-150**	**STOP4-340-150**	**STOP4-425-150**

Prama RF final drills and related stops



h	ø 3.80	ø 4.25	ø 5.00
8.50	**SH-FK380-085**	**SH-FK425-085**	**SH-FK500-085**
10.00	**SH-FK380-100**	**SH-FK425-100**	**SH-FK500-100**
11.50	**SH-FK380-115**	**SH-FK425-115**	**SH-FK500-115**
13.00	**SH-FK380-130**	**SH-FK425-130**	**SH-FK500-130**
15.00	**SH-FK380-150**	**SH-FK425-150**	**SH-FK500-150**



ø 3.80	ø 4.25	ø 5.00
SH-STOP4-FK380	**SH-STOP4-FK425**	**SH-STOP4-FK500**

Countersink drills



ø 3.80	ø 4.25	ø 5.00
A-FCC-380	**A-FCC-425**	**A-FCC-500**

Easy Insert Drivers for implant insertion



EASYC4-EX230-CA short



EASYL4-EX230-CA long



EASY4-EX230-EX with hexagonal connector

Extensions and adapters



BPM-15	extension for instruments with hexagonal connector



PROF-CAL3	extension for surgical drills



AVV-CA-DG-EX	adapter for right angle shank instruments




B-AVV-CA3	mechanical adapter for instruments with hexagonal connector



AVV3-MAN-DG	digital knob

Prama X-ray templates



LA-L100	real dimensions
LA-L120	dimensions increased by 20%
LA-L130	dimensions increased by 30%

Prama RF X-ray templates



LS-L100	real dimensions
LS-L120	dimensions increased by 20%
LS-L130	dimensions increased by 30%

For further information please consult the prosthetic manual MP-IMP-PRAMA-E.

Optional sets

L-INTEGRA-060 set

L-INTEGRA-060 set includes all the necessary components for the insertion of Prama and Prama RF h 6.00 mm. Slots for these instruments are provided in the surgical kit..

Initial drills for Prama h 6.00 mm and related stop		
	**FPS-200**	pilot drill
	**STOPS4-200-060**	stop for pilot drill
	**FGS-200/300**	intermediate drill 200/300
	**FGS-340/425**	intermediate drill 340/425

Instruments included in the ZSHORTY-INT kit.

Final drills for Prama h 6.00 mm and related stops			
	ø 3.80 🟢	ø 4.25 🔵	ø 5.00 🩷
	**FFS-300**	**FFS-340**	**FFS-425**
	**STOPS4-300-060**	**STOPS4-340-060**	**STOPS4-425-060**

Instruments included in the ZSHORTY-INT kit.

Final drills for Prama RF h 6.00 mm and related stops		
	ø 4.25 🔵	ø 5.00 🩷
	**SE-FK410**	**SE-FK500**
	**SE-STOP-FK410-060**	**SE-STOP -FK500-060**

Pins for neck positioning of Prama h 6.00 mm	
	**L-PP-060**

Reply for Prama RF h 6.00 mm	
	**SE-410-060-RP**
	SE-500-060-RP

L-INTEGRA set

L-INTEGRA set includes all the necessary components for the insertion of Prama Slim ø 3.30 mm, the countersink drills and the pins for neck positioning. These components are included in the ultimate version of the Prama surgical kit, but not in the previous version of the surgical kit, thus they are still available for those who have the previous version.

Pins for neck positioning		
	h	
	8.50	**L-PP-085**
	10.00	**L-PP-100**
	11.50	**L-PP-115**
	13.00	**L-PP-130**
	15.00	**L-PP-150**

Countersink drills			
	ø 3.80 🟢	ø 4.25 🔵	ø 5.00 🩷
	A-FCC-380	**A-FCC-425**	**A-FCC-500**

Instruments for Prama Slim implants		
	**FFT3-280-LXS**	final drill
	h	
	8.50	**STOP4-280-085**
	10.00	**STOP4-280-100**
	11.50	**STOP4-280-115**
	13.00	**STOP4-280-130**
	15.00	**STOP4-280-150**
	**A-MS-330**	bone tap with hexagonal connector

Optional instruments (not included in Prama surgical kit)

Intermediate drills

	**FG-200/280XS**	intermediate drill 200/280
	**FG-330/425XS**	intermediate drill 330/425

Drills for distal sectors



ø		ø	
2.00	**FPT5-200-LXS**	3.30	**FFT5-330-LXS**
2.80	**FFT5-280-LXS**	3.40	**FFT5-340-LXS**
2.90	**FFT5-290-LXS**	3.60	**FFT5-360-LXS**
3.00	**FFT5-300-LXS**	4.25	**FFT5-425-LXS**
3.20	**FFT5-320-LXS**	4.55	**FFT5-445-LXS**

Prama bone taps



with right angle shank

ø 3.30	ø 3.80	ø 4.25	ø 5.00
A-MS-330-CA	**A-MS-380-CA**	**A-MS-425-CA**	**A-MS-500-CA**



with hexagonal connector, short


ø 3.30	ø 3.80	ø 4.25	ø 5.00
A-MSC-330	**A-MSC-380**	**A-MSC-425**	**A-MSC-500**

Prama RF bone taps with hexagonal connector



ø 3.80	ø 4.25	ø 5.00
SH-MS-380	**SH-MS-425**	**SH-MS-500**

Optional prosthetic screwdrivers

	**HSMXL-20-EX**	with hexagonal connector, extra long
	**HSMXS-20-CA**	with right angle shank
	**L-HSM-EX**	Full Head, with hexagonal connector, short
	L-HSML-EX	Full Head, with hexagonal connector, long
	L-HSMXL-EX	Full Head, with hexagonal connector, extra long
	**L-HSM-CA**	Full Head, with right angle shank
	**AVV-CAP-TIT-1**	instrument for maintaining of the titanium cap for ball attachments
	**BASCC-EX**	for ball attachments, with hexagonal connector
	**AVV2-ABUT**	for straight P.A.D. abutments
	**8926-SW**	for Locator abutment, short
	**8927-SW**	for Locator abutment, long
	DSPDCLH-24	for Dynamic Abutment, length 24 mm
	**DSPDCLH-32**	for Dynamic Abutment, length 32 mm

Mounter

	**MOU-EX230**	mounter
	**CM2**	mounter stop key

Prama osteotomes



flat tip	
E-OS-020-PP	tip ø 0.20
E-OS-090-PP	tip ø 0.90
concave tip	
E-OS-160-PC	tip ø 1.60
E-OS-200-PC	tip ø 2.00
E-OS-240-PC	tip ø 2.40

Prama RF osteotomes



	flat tip	round tip
h	ø 3.80	
8.50 - 10.00	**SH-OS-380-100-PP**	**SH-OS-380-100-PR**
11.50	**SH-OS-380-115-PP**	**SH-OS-380-115-PR**
13.00	**SH-OS-380-130-PP**	**SH-OS-380-130-PR**
15.00	**SH-OS-380-150-PP**	**SH-OS-380-150-PR**
	ø 4.25	
8.50 - 10.00	**SH-OS-425-100-PP**	**SH-OS-425-100-PR**
11.50	**SH-OS-425-115-PP**	**SH-OS-425-115-PR**
13.00	**SH-OS-425-130-PP**	**SH-OS-425-130-PR**
15.00	**SH-OS-425-150-PP**	**SH-OS-425-150-PR**
	ø 5.00	
8.50 - 10.00	**SH-OS-500-100-PP**	**SH-OS-500-100-PR**
11.50	**SH-OS-500-115-PP**	**SH-OS-500-115-PR**
13.00	**SH-OS-500-130-PP**	**SH-OS-500-130-PR**
15.00	**SH-OS-500-150-PP**	**SH-OS-500-150-PR**

Extraction tools for intraoral removing of implants

	**BC-EX230**	short
	**BL-EX230**	long

Depth gauge



PROF3

P.A.D. carrier



PAD-CAR

Torque wrench with control lever



TWL

Parallelism pins

	**PPTS-2-28**	with depth laser marking, small
	**PP-2/28**	without depth laser marking, standard
	**PPS-2/3**	for short implants

Spare O-ring



ORING180-088

For further information please consult the prosthetic manual MP-IMP-PRAMA-E.

Tissues conditioning

Healing abutments laying on the implant platform

Straight healing abutments



transgingival height	unique platform
2.00	**A-TMG-330-2**
3.00	**A-TMG-330-3**
5.00	**A-TMG-330-5**
7.00	**A-TMG-330-7**

Healing abutments with coronal enlargment



transgingival height	unique platform
2.00	**A-TMTCR-330-2**
3.00	**A-TMTCR-330-3**
5.00	**A-TMTCR-330-5**

Healing abutments with wide profile



transgingival height	unique platform
2.00	**A-TMGB-330-2**
3.00	**A-TMGB-330-3**
4.00	**A-TMGB-330-4**

Temporary posts

Temporary posts

	unique platform	
	**A-PPF-330-EX**	in REEF resin, repositionable
	A-PPF-330	in REEF resin, rotating
	**A-MPSC-330**	Simple in PEEK, with Gr. 5 titanium base, repositionable
	**A-MPSCI-330-EX**	Simple in Gr. 5 titanium base, repositionable
	A-MPSCI-330	Simple in Gr. 5 titanium base, rotating
	**A-CTI-330**	in titanium with castable sleeve, rotating
	A-CCI-S	spare castable sleeve for luting technique

VM2-180 fixation screw included.

Prama IN healing abutments (closing on the implant neck)

Prama IN healing abutments in titanium



transgingival height	closure on the neck	ø 3.80	ø 4.25	ø 5.00
2.50	0.50*	**L-TMG-380-25-05**	**L-TMG-425-25-05**	**L-TMG-500-25-05**
3.50	0.50*	**L-TMG-380-35-05**	**L-TMG-425-35-05**	**L-TMG-500-35-05**
2.50	1.50**	**L-TMG-380-25-15**	**L-TMG-425-25-15**	**L-TMG-500-25-15**
3.50	1.50**	**L-TMG-380-35-15**	**L-TMG-425-35-15**	**L-TMG-500-35-15**

Prama IN healing abutments in REEF resin



transgingival height	closure on the neck	ø 3.80	ø 4.25	ø 5.00
6.00	0.50*	**L-TMGPF-380-05**	**L-TMGPF-425-05**	**L-TMGPF-500-05**
6.00	1.50**	**L-TMGPF-380-15**	**L-TMGPF-425-15**	**L-TMGPF-500-15**

Healing abutments in REEF resin: VM2-180 fixation screw included.
**Only for Regular and Long Neck. | **Only for Regular Neck.*

Impression and model phase

Pick-up transfers



ø 3.30	ø 3.80	ø 4.25	ø 5.00
A-TRA-330	**L-TRA-380**	**L-TRA-425**	**L-TRA-500**

VTRA2-180-15 screw included.

Pull-up transfers



unique platform
A-TRAP-330

Analogs

Short Neck



ø 3.30	ø 3.80	ø 4.25
A-ANA-330	**LS-ANA-380**	**LS-ANA-425**

Regular Neck



ø 3.30	ø 3.80	ø 4.25	ø 5.00
A-ANA-330	**L-ANA-380**	**L-ANA-425**	**L-ANA-500**

Long Neck



ø 3.80	ø 4.25	ø 5.00
LL-ANA-380	**LL-ANA-425**	**LL-ANA-500**

Posts

Premade posts

Straight premade posts			
	h	unique platform	
	6.00	**L-MD-340-6**	straight, repositionable
	8.00	**L-MD-340-8**	
	6.00	**L-MD-340-6-ROT**	conical, rotating
	8.00	**L-MD-340-8-ROT**	

L-VMS-180 fixation screw included.

Angled premade posts			
	transgingival height	unique platform	
	0.00	**L-MA15-340***	angled at 15°, repositionable
	1.00	**L-MA15-340-1**	
	2.00	**L-MA15-340-2**	
	1.00	**L-MA25-340-1**	angled at 25°, repositionable
	2.00	**L-MA25-340-2**	

*L-VMS-180 fixation screw included. | *L-MA15-340 posts: L-VM-180 fixation screw included.*

Preparable posts

Standard preparable posts (laying on the implant platform)		
	unique platform	
	**A-MF-330**	straight emergence
	**A-MFS-330**	wide emergence

VM2-180 fixation screw included.

Prama IN preparable posts (closing on the implant neck)				
	closure on the neck	ø 3.80 🟢	ø 4.25 🔵	ø 5.00 🔴
	0.50*	**L-MF-380-05**	**L-MF-425-05**	**L-MF-500-05**
	1.50**	**L-MF-380-15**	**L-MF-425-15**	**L-MF-500-15**

L-VMS-180 fixation screw included.
** Only for Regular and Long Neck | ** Only for Regular Neck.*

Posts for luting technique

Posts in titanium and in cobalt chrome - rotating		
	unique platform	
	**L-CT-340-ROT**	Gr. 5 titanium posts, h 14.00 mm, straight emergence
	L-CCRCO-340-ROT	Cobalt chrome titanium posts, h 14.00 mm, straight emergence
	**L-CTR-340-ROT**	Gr. 5 titanium posts, h 14.00 mm, anatomical emergence
	L-CCRCOR-340-ROT	Cobalt chrome posts, h 14.00 mm, anatomical emergence

L-VMS-180 fixation screw included.

Castable posts with a metal base and entirely castable posts

Castable abutments with a metal base		
	unique platform	
	**A-UCR-330**	gold alloy base, rotating
	A-UCRCO-330	cobalt chrome base, rotating
	**A-UCRCO-330-EX**	cobalt chrome base, repositionable
	A-UCTR-330-EX	titanium base, repositionable
	A-UCR-330-EX	gold alloy base, repositionable
	**A-CCUCR-330**	spare castable PMMA sleeve, fixation screw not included

VM2-180 fixation screw included.

Dynamic Abutment for angled screw hole		
	unique platform	
	**PD3PKH330/CC**	cobalt chrome base, repositionable
	PD3PKR330/CC	cobalt chrome base, rotating
	**PD3PKH330/P**	entirely castable, repositionable
	PD3PKR330/P	entirely castable, rotating

TPDH18L66 fixation screw not included.

Entirely castable posts		
	unique platform	
	**A-CC-330-EX**	repositionable
	A-CC-330	rotating

VM2-180 fixation screw included.

Prama IN castable posts with a metal base (closing on the implant neck)					
	closure on the neck	ø 3.80 🟢	ø 4.25 🔵	ø 5.00 🔴	
	0.50*	**L-UCRCO-380-05**	**L-UCRCO-425-05**	**L-UCRCO-500-05**	cobalt chrome base, repositionable
	1.50**	**L-UCRCO-380-15**	**L-UCRCO-425-15**	**L-UCRCO-500-15**	cobalt chrome base, repositionable
	0.50*	**L-UCRCO-380-05-ROT**	**L-UCRCO-425-05-ROT**	**L-UCRCO-500-05-ROT**	cobalt chrome base, rotating
	1.50**	**L-UCRCO-380-15-ROT**	**L-UCRCO-425-15-ROT**	**L-UCRCO-500-15-ROT**	cobalt chrome base, rotating
	0.50*	**L-UCR-380-05-ROT**	**L-UCR-425-05-ROT**	**L-UCR-500-05-ROT**	gold alloy base, rotating
		**A-CCUCR-380**	**A-CCUCR-425**	**A-CCUCR-500**	spare castable PMMA sleeve, fixation screw not included

** Only for Regular and Long Neck | ** Only for Regular Neck.*
L-VMS-180 fixation screw included.

Abutment

P.A.D. abutment (Disparallel Screwed Prosthesis)

Straight P.A.D. abutments



transgingival height	unique platform
1.50	**A-PAD-AD330-15**
3.00	**A-PAD-AD330-30**
4.00	**A-PAD-AD330-40**

Angled P.A.D. abutments





transgingival height	unique platform	
3.00	**A-PAD-AA330-173**	angled at 17°
5.00	**A-PAD-AA330-175**	
3.00	**A-PAD-AA330-303**	angled at 30°
5.00	**A-PAD-AA330-305**	

PAD-VM-180 fixation screw included.

P.A.D. components for overstructures

Protection caps





unique platform	
PAD-CG	in Gr. 5 titanium
PAD-CGP	in PEEK

PAD-VP-140 fixation screw included.

Impression and model phase







unique platform	
PAD-CAP	POM cap for impression taking, rotating
PAD-CAP-EX	POM cap for impression taking, repositionable
PAD-TRA	Gr. 5 titanium Pick-up transfer, rotating
PAD-TRA-EX	Gr. 5 titanium Pick-up transfer, repositionable
PAD-ANA	Gr. 5 titanium analog for P.A.D. abutment

Pick-up transfer:PAD-VTRAL-140 fixation screw included.

P.A.D. prostheses for D.P.F. technique











unique platform	
PAD-CT-LV	spare titanium sleeve
PAD-CC-LV	spare castable centring device
PAD-TR-LV	spare anti-escape plug
PAD-ORING-LV	spare protection O-ring
PAD-LV	complete pack of the prosthetic components for the “D.P.F.” technique
BARC	castable bar, l 5 cm, ø 2.20 mm

Titanium sleeves: PAD-VP-140 fixation screw included.

P.A.D. components for prosthesis relining and luting technique









unique platform	
PAD-CP	PEEK sleeve, rotating
PAD-CP-EX	PEEK sleeve, repositionable
PAD-CT	Gr. 5 titanium sleeve, rotating
PAD-CT-EX	Gr. 5 titanium sleeve, repositionable
PAD-CCEM	PMMA castable sleeve for Gr. 5 titanium sleeves, fixation screw not included
PAD-CC	PMMA castable sleeve, rotating
PAD-CC-EX	PMMA castable sleeve, repositionable

Sleeves in PEEK, Gr. 5 titanium and PMMA: PAD-VP-140 fixation screw included.

Castable abutments with a metal base





unique platform	
PAD-UC	PMMA castable sleeve, gold alloy base, rotating
PAD-UCRCO	PMMA castable sleeve, cobalt chrome base, rotating

PAD-VP-140 fixation screw included.

PLAIN abutments

PLAIN abutments



transgingival height	unique platform
2.00	**A-PLAIN-ABU330-2**
3.00	**A-PLAIN-ABU330-3**
4.00	**A-PLAIN-ABU330-4**

Components for PLAIN abutments











unique platform	
A-PLAIN-CG330	healing cap
A-PLAIN-CC330	castable sleeve
A-PLAIN-CT330	titanium sleeve
A-PLAIN-TRA-330	Pick-up transfer
A-PLAIN-ANA-330	analog

Castable and titanium sleeves: A-PLAIN-VP200 fixation screw included.
Pick-up transfer: A-PLAIN-VTRA200 fixation screw included.
For further information please consult the prosthetic manual MP-IMP-PRAMA-E.

CAD-CAM prosthetic solutions

Scanbodies for laboratory and intraoral use

Implant scanbody

	unique platform	
	**L-INT-CAMTRA-L**	h 10.00 mm, in titanium

Implant scanbody: VM2-180 fixation screw included.

Abutment scanbody

	unique platform	
	**PAD-INT-CAMTRA-L**	for P.A.D. abutment, in titanium
	**A-INT-CAMTRA-PLAIN-L**	for PLAIN abutment, in titanium

P.A.D. abutment scanbody: PAD-VP-140 fixation screw included.
PLAIN abutment scanbody: A-PLAIN-VP200 fixation screw included.

Intermediate elements for luting

Standard T-Connect (laying on the implant platform)

	cone height	unique platform	
	4.00	**A-BASTZR-S-330-4**	repositionable
	6.00	**A-BASTZR-S-330-6**	
	4.00	**A-BASTZR-M-330-4**	rotating
	6.00	**A-BASTZR-M-330-6**	

VM2-180 fixation screw included.

T-Connect with different transgingival heights

	transgingival height	cone height	unique platform	
	1.00	4.00	**A-TCR-S-330-4-1**	repositionable
	2.00	4.00	**A-TCR-S-330-4-2**	
	3.00	4.00	**A-TCR-S-330-4-3**	
	1.00	4.00	**A-TCR-M-330-4-1**	rotating
	2.00	4.00	**A-TCR-M-330-4-2**	
	3.00	4.00	**A-TCR-M-330-4-3**	

VM2-180 fixation screw included.

Prama IN T-Connect supports (closing on the implant neck at 0.50 mm)*

	cone height	ø 3.80 🟢	ø 4.25 🔵	ø 5.00 🔴	
	4.00	**L-BASTZR-S-380-4-05**	**L-BASTZR-S-425-4-05**	**L-BASTZR-S-500-4-05**	repositionable
	6.00	**L-BASTZR-S-380-6-05**	**L-BASTZR-S-425-6-05**	**L-BASTZR-S-500-6-05**	
	4.00	**L-BASTZR-M-380-4-05**	**L-BASTZR-M-425-4-05**	**L-BASTZR-M-500-4-05**	rotating
	6.00	**L-BASTZR-M-380-6-05**	**L-BASTZR-M-425-6-05**	**L-BASTZR-M-500-6-05**	

**Only for Regular and Long Neck. | L-VMS-180 fixation screw included.*

Interfase Dynamica with support base and palatal screw hole

	cone height	unique platform	
	4.00	**IND3PKH330/TIA**	repositionable
	4.00	**IND3PKR330/TIA**	rotating

TPDH18L66 fixation screw not included.

Interfase without support base and with palatal screw hole

	cone height	unique platform	
	6.00	**L-MDT-340-6**	repositionable
	6.00	**L-MDT-340-6-ROT**	rotating

TPDH18L66 fixation screw not included.

Sleeves for wax-up modeling

	cone height	ø 3.30 🔵	ø 3.80 🟢	ø 4.25 🔵	ø 5.00 🔴
	4.00	**A-CCBAS-330-4**	**A-CCBAS-330-4**	**A-CCBAS-425-4**	**A-CCBAS-500-4**
	6.00	**A-CCBAS-330-6**	**A-CCBAS-330-6**	**A-CCBAS-425-6**	**A-CCBAS-500-6**

Titanium sleeve for P.A.D. abutment

	**PAD-BASTZR**

PAD-VP-140 fixation screw included.

Chairside T-Connect supports with cementing cone h 5.00 mm

	transgingival height	unique platform
	0.50	**A-TCS-330S-C**
	1.50	**A-TCS-330S-C-15**
	2.50	**A-TCS-330S-C-25**

VM2-180 fixation screw included.

Chairside Prama In T-Connect supports (closing on the implant neck at 0.50 mm)*

	ø 3.80 🟢	ø 4.25 🔵	ø 5.00 🔴
	**L-TCS-380S-05**	**L-TCS-425L-05**	**L-TCS-500L-05**

**Only for Regular and Long Neck. | L-VMS-180 fixation screw included.*

Premilled posts

Titanium Prama IN premilled posts (closing on the implant neck at 0.50 mm)*

	ø 3.80 🟢	ø 4.25 🔵	ø 5.00 🔴
	**L-CAMT-38005-12-SW**	**L-CAMT-42505-12-SW**	**L-CAMT-50005-12-SW**

*VM2-180 fixation screw included. | *Only for Regular and Long Neck.*

Titanium premilled post (laying on the implant platform)

	unique platform
	**A-CAMT-330-12-SW**

VM2-180 fixation screw included.

For further information please consult the prosthetic manual MP-IMP-PRAMA-E.

Conico conometric prosthesis

Straight Conico posts		
	transgingival height	unique platform
	0.50	**A-MD-TS-330-05**
	1.00	**A-MD-TS-330-10**
	2.00	**A-MD-TS-330-20**
	3.00	**A-MD-TS-330-30**

VM2-180 fixation screw included.

Gr. 2 titanium bars, for prostheses on conometric posts*		
	DW-BARRA1.2	circular profile l 150 mm, ø 1.20 mm
	DW-BARRA1.5	circular profile l 150 mm, ø 1.50 mm
	DW-BARRA1.8	circular profile l 150 mm, ø 1.80 mm
	DW-BARRA1x3	rectangular profile l 100 mm, 3x1 mm
	DW-BARRA2x4	rectangular profile l 100 mm, 4x2 mm

** 5 pcs packaging.*

Angled Conico posts			
	transgingival height	unique platform	
	2.00	**A-MA07-TS-330-2**	angled at 7.5°
	2.00	**A-MA15-TS-330-2**	angled at 15°
	2.00	**A-MA20-TS-330-2**	angled at 22.5°

VM2-180 fixation screw included.

Conometric components		
	**CAP2-TS-DEF**	cap for partial or full-arch fixed prosthesis
	**CAP2-TS-IND**	non-rotating cap for single element
	**CAP2-TS-REM**	cap for the removal
	**CAP2-TS-IMP**	PEEK transfer cap for direct impression on Conico abutment
	**ANA2-MD-TS**	Conico post analog
	**CAP2-TS-PIN**	cap with pin, to be used for intraoral welding

Accessories for Conico system

	**PAR-PP**	parallelometer		**PAR-INS-AS**	parallelometer insert

Remember to always purchase the parallelometer with the insert relating to the implant connection on which it is to be used.

Clinical and laboratory prosthetic screws

Standard fixation screws		
	**VM2-180**	
	VM2-180-10	10 pcs packaging

Conical support fixation screws		
	**L-VMS-180**	
	L-VMS-180-10	10 pcs packaging

Full Head technology fixation screws		
	**L-VM-180**	
	L-VM-180-10	10 pcs packaging

Fixation screws for Dynamic Abutment and Interfase Dynamica	
	**TPDH18L66**

Screws for Pick-up transfer		
	**VTRA2-180-15**	h 15.00
	**VTRA2-180-20**	h 20.00
	**VTRA2-180-MAN**	for manual screwing

Recommended torque for tightening prosthetic components	
8-10 Ncm	surgical cover screws, healing abutments
20-25 Ncm	prosthetic screws
25-30 Ncm	direct screwing on implants prosthetic components
8-10 Ncm	transfer screws

Screws for P.A.D. abutments		
	**PAD-VM-180**	
	PAD-VM-180-10	10 pcs packaging
	**PAD-VTRAL-140**	for Pick-up transfer, long
	**PAD-VTRA-140**	for Pick-up transfer, standard
	**PAD-VTRAL-140-MAN**	for Pick-up transfer, manual screwing
	**PAD-VP-140**	spare screw
	PAD-VP-140-10	10 pcs packaging

Screws for PLAIN abutments		
	**A-PLAIN-VP200**	for PLAIN abutments
	A-PLAIN-VP200-10	10 pcs packaging
	**A-PLAIN-VTRA200**	for PLAIN abutment transfers

Scews for CAD-CAM custom made prosthesis		
	**A-CAMTVABU180**	for zirconium oxide single crowns and ECHO overstructures
	**PAD-VCAM140**	for zirconium oxide ECHO overstructures on P.A.D. abutments
	**CAMPRON205-10**	spare ring washers for the head of the screw

Warning
Fixation screws are always included in the packaging of each prosthetic component, unless otherwise indicated.
Please use laboratory screws for laboratory phases and keep the screw included with the prosthetic component for the final fastening in the patient's mouth.

Solutions for overdenture

Locator and Locator R-Tx abutments

Abutments



transgingival height	Locator	Locator R-Tx
1.00	**1670**	**31110-01**
2.00	**1671**	**31110-02**
3.00	**1672**	**31110-03**
4.00	-	**31110-04**
5.00	**1674**	**31110-05**
6.00	-	**31110-06**



Accessories and spare parts for Locator abutments

	unique platform	
	**8530**	4 aluminium analogs
	**8505**	4 aluminium transfers, 4 retainers cod. 8515 included
	**8517**	4 parallelism pins in nylon
	**9530**	steel plate AISI 316 L

Locator Core Tool

	**8393**	complete Locator Core Tool
	**8397**	steel push rod for inserting the retainers in the caps
	**8390**	steel driver for abutment screwing/ unscrewing
	8394	retention jacket for the driver, to transfer the Locator abutments into the oral cavity

Accessories for Locator R-Tx abutment

	unique platform	
	**30013**	pack of 4 caps and 4 black retainers for impressione taking
	**30001**	grey, without retention, pack of 4 pcs
	**30002**	blue, low retention, pack of 4 pcs
	**30003**	pink, medium retention, pack of 4 pcs
	**30004**	white, high retention, pack of 4 pcs
	**30012**	black, for the impression taking, pack of 4 pcs
	**8514**	spacer rings for the prosthetis relining phase, pack of 20 pcs
	30014	aluminium analogs for abutments ø 3.35 mm pack of 4 pcs
	30015	aluminium analogs for abutments ø 4.00 mm, pack of 4 pcs
	30016	aluminium analogs for abutments ø 5.00 mm, pack of 4 pcs
	**30017**	transfers with low retention, pack of 4 pcs
	30018	spacer rings for relining phase, pack of 4 pcs
	**30021**	Locator Core Tool for Locator R-Tx
	08008	RA screwdriver for Locator R-Tx, hex 1.25 mm
	**08007**	digital screwdriver for Locator R-Tx, hex 1.25 mm

Accessories and spare parts for Locator abutments

	unique platform	
	**8515**	low retention capacity for impression taking, black, in polyethylene, 4 pcs packaging
	**8524**	retention 5 lb = 2268 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	**8527**	retention 3 lb = 1361 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	**8529**	retention 1,5 lb= 680 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	**8547**	retention 4 lb= 1814 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	**8915**	retention 2 lb= 907 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	**8548**	retention 1 lb= 450 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	8514	spacer rings in silicon rubber, 20 pcs packaging
	**8519-2**	kit: 2 Gr. 5 titanium caps, 2 spacer rings and 2 nylon retainers for each of the 4 different retention capacities
	**8540-2**	kit: 2 Gr. 5 titanium caps, 2 spacer and 2 nylon retainers for each of the 4 different retention capacities meant for disparallalelisms
	**8550-2**	kit: 2 Gr. 5 steel caps, 2 spacer rings and 2 nylon retainers for each of the 4 different retention capacities

Ball attachments and accessories

Ball attachments



transgingival height	unique platform
1.00	**A-AS-330-1**
2.00	**A-AS-330-2**
4.00	**A-AS-330-4**

Analog for ball attachments



unique platform
ANAS

Caps for ball attachments


	unique platform	
	**CAP-TFL-1**	polyamide cap
	**CONT-CAP-TFL-1**	steel container for polyamide cap
	CAP-TIT-1	Gr. 5 Ti cap with cap in two parts, Ti retention spring and plastic
	AN-CAP-TIT-1	spare plastic ring
	MOL1-CAP-TIT-1	spare steel spring, average hardness
	MOL2-CAP-TIT-1	spare steel spring, soft
	CAP-1	gold alloy cap, with 3 plastic rings and a tin spacer for the lab

Components for bars

Overdenture on bars

	BARC-CAV-TIT	castable bar, l 5.00 cm, h 3.00 mm, thickness 2.20 mm, ovoid-shaped profile with spacer
	CAV-TIT	divisible bar attachment in titanium for oval bars, h 3.00 mm, thickness 2.20 mm
	BARC	castable bar, l 5.00 cm, ø 2.20 mm
	CAV-375	bar attachment in gold alloy for round bars ø 2.20 mm

For further information please consult the prosthetic manual MP-IMP-PRAMA-E.



FUTURE
DENTAL

FUTURE DENTAL
България, гр. София
кв. Кръстова вада
ул. „ген. Иван Колев“ 14
тел.: 088 3333 730