

PREMIUM ONE

Product overview



FUTURE DENTAL

Premium One implants



h	ø 3.30	ø 3.80	ø 4.25	ø 5.00
7.00	-	-	**AS-ZT-425-070**	**AS-ZT-500-070**
8.50	**A-ZT-330-085**	**A-ZT-380-085**	**AS-ZT-425-085**	**AS-ZT-500-085**
10.00	**A-ZT-330-100**	**A-ZT-380-100**	**AS-ZT-425-100**	**AS-ZT-500-100**
11.50	**A-ZT-330-115**	**A-ZT-380-115**	**AS-ZT-425-115**	**AS-ZT-500-115**
13.00	**A-ZT-330-130**	**A-ZT-380-130**	**AS-ZT-425-130**	**AS-ZT-500-130**
15.00	**A-ZT-330-150**	**A-ZT-380-150**	**AS-ZT-425-150**	**AS-ZT-500-150**
18.00	-	**A-ZT-380-180**	**AS-ZT-425-180**	-

Surgical cover screw				
	**A-VT-330**	**A-VT-380**	**SH-VT-425-BL**	**SH-VT-500-VI**

Surgical kits

Premium One surgical kit



ZPREMIUM-ONE-INT	grommetless surgical kit including all the instruments for Premium One implants
AS-TRAY-INT	Radel instrument grommetless tray for Premium One implants


For further information on the surgical kit and the surgical procedures, please consult the surgical manual MC-IMP-PREMIUM-ONE-E.

Drilling kit for short Premium One implants



ZSHORTY-INT	drilling kit for short Premium One implants
SHORTY-TRAY-INT	Radel instrument grommetless tray for short Premium One implants

Premium One and Shelta surgical kit



ZPRESH-INT	grommetless surgical kit including all the instruments for Premium One and Shelta implants
ASH-TRAY-INT	Radel instrument grommetless tray for Premium One and Shelta implants

For further information on the surgical kit and the surgical procedures, please consult the surgical manual MC-IMP-

Premium Kohno One surgical kit



ZPREKOH-ONE-INT	grommetless surgical kit including all the instruments for Premium One and Kohno One implants
AK-TRAY-INT	Radel instrument grommetless tray for Premium One and Kohno One implants

Premium One guided surgery surgical kit



ZGS-PRO-A-INT	grommetless surgical kit including all the instruments for guided surgery insertion of Premium One implants
GSPROA-TRAY-INT	Radel instrument grommetless tray for guided insertion of Premium One implants

For further information on the surgical kit and the surgical procedures, please consult the surgical manual MC-UD-EPLANPROA-E.

Instruments included in Premium One surgical kit

Precision and intermediate drills

	**FS-230**	precision drill
	**FG-200/280XS**	intermediate drill 200/280
	**FG-330/425XS**	intermediate drill 330/425

Pilot drill, intermediate drill and related stops

	FPT3-200-LXS	pilot drill
	**FFT3-250-LXS**	intermediate cylindrical drill

	h	
	7.00	**STOP4-200-070**
	8.50	**STOP4-200-085**
	10.00	**STOP4-200-100**
	11.50	**STOP4-200-115**
	13.00	**STOP4-200-130**
	15.00	**STOP4-200-150**

Bone taps with hexagonal connector

	ø 3.30	ø 3.80	ø 4.25	ø 5.00
	**A-MS-330**	**A-MS-380**	**A-MS-425**	**A-MS-500**

Screwdrivers

	**HSMXS-20-DG**	digital, extra short
	HSM-20-DG	digital, short
	HSML-20-DG	digital, long
	**HSM-20-EX**	with hexagonal connector
	HSML-20-EX	with hexagonal connector, long
	**HSM-20-CA**	with right angle shank

Dynamometric ratchet



CRI5-KIT

Easy Insert Drivers for implant insertion

	**EASYC4-EX230-CA**	short
	**EASYL4-EX230-CA**	long
	**EASY4-EX230-EX**	with hexagonal connector

Final drills and related stops

	ø 3.30	ø 3.80	ø 4.25	ø 5.00
	**FFT3-280-LXS**	**FFT3-300-LXS**	**FFT3-340-LXS**	**FFT3-425-LXS**

	h	ø 3.30	ø 3.80	ø 4.25	ø 5.00
	7.00	**STOP4-280-070**	**STOP4-300-070**	**STOP4-340-070**	**STOP4-425-070**
	8.50	**STOP4-280-085**	**STOP4-300-085**	**STOP4-340-085**	**STOP4-425-085**
	10.00	**STOP4-280-100**	**STOP4-300-100**	**STOP4-340-100**	**STOP4-425-100**
	11.50	**STOP4-280-115**	**STOP4-300-115**	**STOP4-340-115**	**STOP4-425-115**
	13.00	**STOP4-280-130**	**STOP4-300-130**	**STOP4-340-130**	**STOP4-425-130**
	15.00	**STOP4-280-150**	**STOP4-300-150**	**STOP4-340-150**	**STOP4-425-150**

Extensions and adapters

	**BPM-15**	extensions for instruments with hexagonal connector
	**PROF-CAL3**	extension for surgical drills
	**B-AVV-CA3**	mechanical adapter for instruments with hexagonal connector
	**AVV3-MAN-DG**	digital knob

Parallelism pins with depth laser marking



PPTL-2-28

X-ray templates

	**AS-L100**	real dimensions
	AS-L120	dimensions increased by 20%
	AS-L130	dimensions increased by 30%

Optional instruments (not included in Premium One surgical kit)

Initial drills for short implants

	**FPS-200**	pilot drill
	**STOPS4-200-070**	stop for pilot drill
	**FGS-200/300**	intermediate drill 200/300
	**FGS-340/425**	intermediate drill 340/425
	**FGS-425/540**	intermediate drill 425/540

Final drills for short implants and related stops

ø 3.80	ø 4.25	ø 5.00
 **FFS-300**	**FFS-340**	**FFS-425**
 **STOPS4-300-070**	**STOPS4-340-070**	**STOPS4-425-070**


Short bone taps with hexagonal connector

ø 3.30	ø 3.80	ø 4.25	ø 5.00
 **A-MSC-330**	**A-MSC-380**	**A-MSC-425**	**A-MSC-500**



Bone taps with right angle shank

ø 3.30	ø 3.80	ø 4.25	ø 5.00
 **A-MS-330-CA**	**A-MS-380-CA**	**A-MS-425-CA**	**A-MS-500-CA**

Screwdrivers

	**HSMXL-20-EX**	with hexagonal connector, extra long
	**HSMXS-20-CA**	with right angle shank, extra short
	**AVV-CAP-TIT-1**	instrument for maintaining of the titanium cap for ball attachments
	**BASCC-EX**	for ball attachments, with hexagonal connector
	**AVV2-ABUT**	for straight P.A.D. abutments
	**8926-SW**	for Locator abutment, short
	**8927-SW**	for Locator abutment, long
	DSPDCLH-24	for Dynamic Abutment and Interfase Dynamica, length 24 mm
	DSPDCLH-32	for Dynamic Abutment and Interfase Dynamica, length 32 mm

Parallelism pins

	**PPTS-2-28**	with depth laser marking, small
	**PP-2/28**	without depth laser marking, standard
	**PPS-2/3**	for short implants

Adapter

	**AVV-CA-DG-EX**	adapter for right angle shank instruments

Countersink drills

ø 3.80	ø 4.25	ø 5.00
 **A-FCC-380**	**A-FCC-425**	**A-FCC-500**

Drills for distal sectors



ø		ø	
2.00	**FPT5-200-LXS**	3.30	**FFT5-330-LXS**
2.80	**FFT5-280-LXS**	3.40	**FFT5-340-LXS**
2.90	**FFT5-290-LXS**	3.60	**FFT5-360-LXS**
3.00	**FFT5-300-LXS**	4.25	**FFT5-425-LXS**
3.20	**FFT5-320-LXS**	4.45	**FFT5-445-LXS**

Bone profilers

transgingival h 3.00, wide



ø 3.30	ø 3.80	ø 4.25	ø 5.00
A-PAD-PS330-L	**A-PAD-PS380-L**	**A-PAD-PS425-L**	**A-PAD-PS500-L**

transgingival h 5.00, narrow

ø 3.30	ø 3.80	ø 4.25	ø 5.00
A-PAD-PS330-S	**A-PAD-PS380-S**	**A-PAD-PS425-S**	**A-PAD-PS500-S**

 **A-PAD-GUI-PS-230** guide cylinder

Osteotomes with flat tip



E-OS-020-PP	ø tip 0.20
E-OS-090-PP	ø tip 0.90
E-OS-160-PP	ø tip 1.60
E-OS-200-PP	ø tip 2.00
E-OS-240-PP	ø tip 2.40

Extraction tools for intraoral removing of implants

	**BC-EX230**	corta
	**BL-EX230**	lunga

Depth gauge

 **PROF3**

Mounter

	**MOU-EX230**	mounter
	**CM2**	mounter stop key

P.A.D. carrier

 **PAD-CAR**

Torque wrench with control lever

 **TWL**

Spare O-ring

 **ORING180-088**

For further information please consult the surgical manual MC-IMP-PREMIUM-ONE-E.

Implant platforms and prosthetic combinations

The connection platforms of Premium One implants are characterized by the presence of the same hexagon and external collar, which allows the choice of the size of the abutment.
According to the chosen prosthetic protocol, a Platform Switching can be made or a matching size post can be used.

ø implant	3.30 ●	3.80 ●	4.25 ●	5.00 ●
ø compatible prosthetic component	3.30 ●	3.30 ●	3.30 ●	3.30 ●
	-	3.80 ●	3.80 ●	3.80 ●
	-	-	4.25 ●	4.25 ●
	-	-	-	5.00 ●

Please note
Prosthetic components with ø 3.30 mm determine prosthetic Platform Switching with ø 3.80 mm implants. These components should only be used for single crowns in frontal sectors (premolar excluded), and in distal sectors only to support multiple restorations.

Tissues conditioning

Healing abutments

Healing abutments with straight emergence



transgingival h	ø 3.30 ●	ø 3.80 ●
2.00	**A-TMG-330-2**	**A-TMG-380-2**
3.00	**A-TMG-330-3**	**A-TMG-380-3**
5.00	**A-TMG-330-5**	**A-TMG-380-5**
7.00	**A-TMG-330-7**	-

Healing abutments with anatomical emergence



transgingival h	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
2.00	**A-TMGR-330-2**	**A-TMGR-380-2**	**AS-TMGR-425-2**	**AS-TMGR-500-2**
3.00	**A-TMGR-330-3**	**A-TMGR-380-3**	**AS-TMGR-425-3**	**AS-TMGR-500-3**
5.00	**A-TMGR-330-5**	**A-TMGR-380-5**	**AS-TMGR-425-5**	**AS-TMGR-500-5**
7.00	-	**A-TMGR-380-7**	**AS-TMGR-425-7**	**AS-TMGR-500-7**

Healing abutments with coronal enlargement



transgingival h	ø 3.30 ●	ø 3.80 ●
2.00	**A-TMTCR-330-2**	**A-TMTCR-380-2**
3.00	**A-TMTCR-330-3**	**A-TMTCR-380-3**
5.00	**A-TMTCR-330-5**	**A-TMTCR-380-5**

Healing abutments with wide profile



transgingival h	ø 3.30 ●	ø 3.80 ●
2.00	**A-TMGB-330-2**	**A-TMGB-380-2**
3.00	**A-TMGB-330-3**	**A-TMGB-380-3**
4.00	**A-TMGB-330-4**	**A-TMGB-380-4**

Temporary posts

PEEK temporary posts with titanium base					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-MPSC-330	-	-	-	repositionable, straight emergence
	A-MPSCR-330	**A-MPSCR-380**	**AS-MPSCR-425**	**AS-MPSCR-500**	repositionable, anatomical emergence

VM2-180 fixation screw included.

SIMPLE titanium temporary posts					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	**A-MPSA-330-EX**	**A-MPSA-380-EX**	-	-	repositionable, anatomical emergence
	A-MPSA-330	**A-MPSA-380**	-	-	rotating, anatomical emergence
	**A-MPSCI-330-EX**	**A-MPSCI-380-EX**	**AS-MPSCI-425-EX**	**AS-MPSCI-500-EX**	repositionable, straight emergence
	A-MPSCI-330	**A-MPSCI-380**	**AS-MPSCI-425**	**AS-MPSCI-500**	rotating, straight emergence
	**A-MPS-330**	**A-MPS-380**	-	-	rotating, wide emergence

VM2-180 fixation screw included.

Titanium temporary posts with castable sleeve					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	**A-CTI-330**	**A-CTI-380**	-	-	rotating
	A-CCI-S	**A-CCI-S**	-	-	spare castable PMMA sleeve, fixation screw not included

VM2-180 fixation screw included.

Impression and model phase

Pick-up transfers					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	**A-TRA-330**	**A-TRA-380**	**AS-TRA-425**	**AS-TRA-500**	straight emergence
	A-TRAR-330	**A-TRAR-380**	**AS-TRAR-425**	**AS-TRAR-500**	anatomical emergence

VTRA2-180-15 screw included.

Pull-up transfers					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	**A-TRAP-330**	-	-	-	straight emergence
	A-TRARP-330	**A-TRARP-380**	**AS-TRARP-425**	**AS-TRARP-500**	anatomical emergence

Closed-tray transfers					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	**A-TRAS-330**	-	-	-	straight emergence
	A-TRARS-330	**A-TRARS-380**	**AS-TRARS-425**	**AS-TRARS-500**	anatomical emergence

VTRA2-180-10 fixation screw included.

Analogs		
	**A-ANA-330**	ø 3.30
	**A-ANA-380**	ø 3.80
	**SH-ANA-425**	ø 4.25
	**SH-ANA-500**	ø 5.00

For further information please consult the prosthetic manual MP-IMP-PKS-ONE-E.

Posts

Premade posts

Premade straight posts						
	transgingival h	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	1.00	**A-MD-330-1**	**A-MD-380-1**	**AS-MD-425-1**	**AS-MD-500-1**	straight emergence
	2.00	**A-MD-330-2**	**A-MD-380-2**	**AS-MD-425-2**	**AS-MD-500-2**	
	4.00	**A-MD-330-4**	**A-MD-380-4**	**AS-MD-425-4**	**AS-MD-500-4**	
	1.00	**A-MDR-330-1**	**A-MDR-380-1**	**AS-MDR-425-1**	**AS-MDR-500-1**	anatomical emergence
	2.00	**A-MDR-330-2**	**A-MDR-380-2**	**AS-MDR-425-2**	**AS-MDR-500-2**	
	4.00	**A-MDR-330-4**	**A-MDR-380-4**	**AS-MDR-425-4**	**AS-MDR-500-4**	

VM2-180 fixation screw included.

Shoulderless premade straight posts						
	h post	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	15.00	**A-MD-330-15**	**A-MD-380-15**	-	-	repositionable

VM2-180 fixation screw included.

Premade angled posts						
	transgingival h	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	1.75	**A-MA15-330**	-	-	-	angled at 15°, straight emergence
	1.80	**A-MAR15-330**	**A-MAR15-380**	**AS-MAR15-425**	**AS-MAR15-500**	angled at 15°, anatomical emergence
	1.80	-	**A-MAR25-380**	**AS-MAR25-425**	-	angled at 25°, anatomical emergence

VM2-180 fixation screw included.

Preparable posts

Standard preparable posts					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-MF-330	**A-MF-380**	**AS-MF-425**	**AS-MF-500**	straight emergence
	A-MFR-330	**A-MFR-380**	**AS-MFR-425**	**AS-MFR-500**	anatomical emergence

VM2-180 fixation screw included.

Precut preparable posts					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-MFP-330	-	-	-	straight emergence
	A-MFPR-330	**A-MFPR-380**	**AS-MFPR-425**	**AS-MFPR-500**	anatomical emergence

VM2-180 fixation screw included.

SIMPLE preparable posts					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-MFS-330	**A-MFS-380**	**AS-MFS-425**	**AS-MFS-500**	very wide emergence

VM2-180 fixation screw included.

Prosthetic solutions for vertical technique

Componenti for vertical technique		
	unique platform	
	**A-TMG-MEFL-330**	Gr. 5 titanium healing abutments
	**A-PPF-330-EX**	REEF resin posts, repositionable
	A-PPF-330	REEF resin posts, rotating
	**A-MEFL-330**	Gr. 5 titanium preparable posts
	**CAP-MEFL-5**	impression cap

Resin and preparable posts: fixation screw VM2-180 included.

XA prosthetic range

Posts for XA cemented restorations



transgingival h	ø 3.30	ø 3.80	ø 4.25	ø 5.00
1.00	SH-MD-F-330-1	SH-MD-F-380-1	SH-MD-F-425-1	-
2.00	SH-MD-F-330-2	SH-MD-F-380-2	SH-MD-F-425-2	-

L-VMS-180 fixation screw included.

Components for XA cemented restorations

	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	SH-ANA-MD-F-330	SH-ANA-MD-F-380	SH-ANA-MD-F-425	-	analogs
	SH-TRA-MD-F-330	SH-TRA-MD-F-380	SH-TRA-MD-F-425	-	transfers
	SH-CP-MD-F-330	SH-CP-MD-F-380	SH-CP-MD-F-425	-	PEEK caps

Posts for screw retained restorations for XA technique



transgingival h	ø 3.30	ø 3.80
4.50	A-ABU-F-TS-330-4	SH-ABU-F-TS-380-4
5.50	A-ABU-F-TS-330-5	SH-ABU-F-TS-380-5
6.50	A-ABU-F-TS-330-6	SH-ABU-F-TS-380-6

Fixation screw included (see details at page 15).

Scanbodies for laboratory and intraoral use in titanium



unique platform	
XA-INT-CAMTRA-L	per pilastri XA

A-PLAIN-VP200 screw included.

Titanium premilled posts for XA abutment



unique platform
XA-CAMT-ABU-12-SW

A-PLAIN-VP100 screw included.

Components for XA screw retained restorations



unique platform	
SH-ANABU-F-380	analogs
SH-TRABU-F-380	transfers
SH-CCABU-F-380	castable sleeves, repositionable
SH-CCABU-F-380-ROT	castable sleeves, rotating
SH-CG-ABU-F-380	PEEK healing cap
SH-CTABU-F-380	Gr. 5 titanium sleeves, repositionable
SH-CTABU-F-380-ROT	Gr. 5 titanium sleeves, rotating
CAP-TS-DEF	Conoweld cap for luting

Transfer: SH-VTRABU-F-200 screw included. Castable and titanium sleeves: A-PLAIN-VP200 screw included.

Castable posts with a metal base and fully castable posts

Castable posts with a metal base

	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-UCR-330-EX	A-UCR-380-EX	AS-UCR-425-EX	AS-UCR-500-EX	gold alloy base, repositionable
	A-UCTR-330-EX	A-UCTR-380-EX	-	-	Gr. 5 titanium base, repositionable
	A-UCRCO-330-EX	A-UCRCO-380-EX	AS-UCRCO-425-EX	AS-UCRCO-500-EX	cobalt chrome base, repositionable
	A-UCR-330	A-UCR-380	AS-UCR-425	AS-UCR-500	gold alloy base, rotating
	A-UCRCO-330	A-UCRCO-380	AS-UCRCO-425	AS-UCRCO-500	cobalt chrome base, rotating
	A-CCUCR-330	A-CCUCR-380	AS-CCUCR-425	AS-CCUCR-500	spare castable PMMA sleeve, fixation screw not included

VM2-180 fixation screw included.

Dynamic Abutment

	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	PD3PKH330/CC	PD3PKH380/CC	-	-	cobalt chrome base, repositionable
	PD3PKR330/CC	PD3PKR380/CC	-	-	cobalt chrome base, rotating
	PD3PKH330/P	PD3PKH380/P	-	-	fully castable, repositionable
	PD3PKR330/P	PD3PKR380/P	-	-	fully castable, rotating

TPDH18L66 screw not included.

Fully castable posts

	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-CC-330-EX	A-CC-380-EX	-	-	repositionable, straight emergence
	A-CCR-330-EX	A-CCR-380-EX	AS-CCR-425-EX	AS-CCR-500-EX	repositionable, anatomical emergence
	A-CC-330	A-CC-380	-	-	rotating, straight emergence
	A-CCR-330	A-CCR-380	AS-CCR-425	AS-CCR-500	rotating, anatomical emergence

VM2-180 fixation screw included.

For further information please consult the prosthetic manual MP-IMP-PKS-ONE-E.

Abutments

P.A.D. abutments (Disparallel Screwed Prosthesis)

Straight P.A.D. abutments



transgingival h	ø 3.30	ø 3.80	ø 4.25	ø 5.00
1.50	**A-PAD-AD330-15**	**A-PAD-AD380-15**	**AS-PAD-AD425-15**	**AS-PAD-AD500-15**
3.00	**A-PAD-AD330-30**	**A-PAD-AD380-30**	**AS-PAD-AD425-30**	**AS-PAD-AD500-30**
4.00	**A-PAD-AD330-40**	**A-PAD-AD380-40**	**AS-PAD-AD425-40**	**AS-PAD-AD500-40**

Angled P.A.D. abutments



transgingival h	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
3.00	**A-PAD-AA330-173**	**A-PAD-AA380-173**	**AS-PAD-AA425-173**	**AS-PAD-AA500-173**	angled at 17°
5.00	**A-PAD-AA330-175**	**A-PAD-AA380-175**	**AS-PAD-AA425-175**	**AS-PAD-AA500-175**	
3.00	**A-PAD-AA330-303**	**A-PAD-AA380-303**	**AS-PAD-AA425-303**	**AS-PAD-AA500-303**	angled at 30°
5.00	**A-PAD-AA330-305**	**A-PAD-AA380-305**	**AS-PAD-AA425-305**	**AS-PAD-AA500-305**	

PAD-VM-180 fixation screw included.

P.A.D. components for overstructures

Protection caps

	unique platform	
	**PAD-CG**	in Gr. 5 titanium
	**PAD-CGP**	in PEEK

PAD-VP-140 fixation screw included.

Impression and model phase

	unique platform	
	**PAD-CAP**	POM cap for impression taking, rotating
	PAD-CAP-EX	POM cap for impression taking, repositionable
	**PAD-TRA**	Gr. 5 titanium Pick-up transfer, rotating
	PAD-TRA-EX	Gr. 5 titanium Pick-up transfer, repositionable
	**PAD-ANA**	Gr. 5 titanium analog for P.A.D. abutment

Pick-up transfer: PAD-VTRAL-140 fixation screw included.

Castable abutments with a metal base

	unique platform	
	**PAD-UC**	PMMA castable sleeve, gold alloy base, rotating
	**PAD-UCRCO**	PMMA castable sleeve, cobalt chrome base, repositionable

PAD-VP-140 fixation screw included.

Scanbodies for laboratory and intraoral use in titanium

	unique platform	
	**PAD-INT-CAMTRA-L**	for P.A.D. abutment

A-PLAIN-VP200 screw included.

P.A.D. components for prosthesis relining and luting technique

	piattaforma unica	
	**PAD-CP**	PEEK sleeve for relining, rotating
	PAD-CP-EX	PEEK sleeve for relining, repositionable
	**PAD-CT**	Gr. 5 titanium sleeves, rotating
	PAD-CT-EX	Gr. 5 titanium sleeves, repositionable
	**PAD-CCEM**	PMMA castable sleeve for Gr. 5 titanium sleeves, fixation screw not included
	**PAD-CC**	castable sleeves in PMMA, rotating
	PAD-CC-EX	castable sleeves in PMMA, repositionable

PEEK, Gr. 5 titanium and PMMA sleeves: PAD-VP-140 screw included.

Titanium sleeves for P.A.D. abutments

	**PAD-BASTZR**

PAD-VP-140 screw included.

P.A.D. prostheses for D.P.F. technique

	unique platform	
	**PAD-CT-LV**	spare titanium sleeve
	**PAD-CC-LV**	spare castable centring device
	**PAD-TR-LV**	spare anti-escape plug
	**PAD-ORING-LV**	spare protection O-ring
	**PAD-LV**	complete pack of the prosthetic components for the “D.P.F.” technique
	BARC	castable bars, l 5 cm, ø 2.20 mm

Titanium sleeves: PAD-VP-140 fixation screw included.

Posts for luting technique

Posts in titanium and in cobalt chrome - rotating

	unique platform	
	**L-CT-340-ROT**	Gr. 5 titanium posts, h 14.00 mm, straight emergence
	L-CCRCO-340-ROT	Cobalt chrome titanium posts, h 14.00 mm, straight emergence

L-VMS-180 screw included.

Posts in titanium and in cobalt chrome - rotating

	unique platform	
	**L-CTR-340-ROT**	Gr. 5 titanium posts, h 14.00 mm, anatomical emergence
	L-CCRCOR-340-ROT	Cobalt chrome posts, h 14.00 mm, anatomical emergence

L-VMS-180screw included.

PLAIN abutments

PLAIN abutments					
	h transmucosa	ø 3.30	ø 3.80	ø 4.25	ø 5.00
	2.00	**A-PLAIN-ABU330-2**	**A-PLAIN-ABU380-2**	**AS-PLAIN-ABU425-2**	**AS-PLAIN-ABU500-2**
	3.00	**A-PLAIN-ABU330-3**	**A-PLAIN-ABU380-3**	**AS-PLAIN-ABU425-3**	**AS-PLAIN-ABU500-3**
	4.00	**A-PLAIN-ABU330-4**	**A-PLAIN-ABU380-4**	**AS-PLAIN-ABU425-4**	**AS-PLAIN-ABU500-4**

Components for PLAIN abutments					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	**A-PLAIN-CG330**	**A-PLAIN-CG380**	**A-PLAIN-CG425**	**A-PLAIN-CG500**	healing caps
	**A-PLAIN-CC330**	**A-PLAIN-CC380**	**A-PLAIN-CC425**	**A-PLAIN-CC500**	castable sleeves
	**A-PLAIN-CT330**	**A-PLAIN-CT380**	**A-PLAIN-CT425**	**A-PLAIN-CT500**	Gr. 5 titanium sleeves
	**A-PLAIN-TRA-330**	**A-PLAIN-TRA-380**	**A-PLAIN-TRA-425**	**A-PLAIN-TRA-500**	Pick-up transfers
	**A-PLAIN-ANA-330**	**A-PLAIN-ANA-380**	**A-PLAIN-ANA-425**	**A-PLAIN-ANA-500**	analogs

Castable and titanium sleeves: A-PLAIN-VP200 fixation screw included. Pick-up transfers: A-PLAIN-VTRA200 fixation screw included.

Scanbodies for laboratory and intraoral use in titanium		
	unique platform	
	A-INT-CAMTRA-PLAIN-L	for PLAIN abutments

A-PLAIN-VP200 screw included.

Conico conometric prosthesis

Straight Conico posts		
	transgingival h	ø 3.30	ø 3.80
	0.50	**A-MD-TS-330-05**	-
	1.00	**A-MD-TS-330-10**	**A-MD-TS-380-10**
	2.00	**A-MD-TS-330-20**	**A-MD-TS-380-20**
	3.00	**A-MD-TS-330-30**	**A-MD-TS-380-30**

VM2-180 fixation screw included.

Angled Conico posts			
	ø 3.30	ø 3.80	
	**A-MA07-TS-330-2**	**A-MA07-TS-380-2**	angolati a 7.5°
	A-MA15-TS-330-2	**A-MA15-TS-380-2**	angolati a 15°
	A-MA20-TS-330-2	**A-MA20-TS-380-2**	angolati a 22.5°

VM2-180 fixation screw included.

Components for conometric prosthesis		
	**CAP2-TS-DEF**	cap for partial or full-arch fixed prosthesis
	**CAP2-TS-IND**	non-rotating cap for single element
	**CAP2-TS-REM**	cap for the removal
	**CAP2-TS-IMP**	PEEK transfer cap for direct impression on Conico abutment
	**ANA2-MD-TS**	Conico post analog
	**CAP2-TS-PIN**	cap with pin, to be used for intraoral welding

Gr. 2 titanium bars, for prostheses on conometric posts*		
	**DW-BARRA1.2**	circular profile l 150 mm, ø 1.20 mm
	**DW-BARRA1.5**	circular profile l 150 mm, ø 1.50 mm
	**DW-BARRA1.8**	circular profile l 150 mm, ø 1.80 mm
	**DW-BARRA1x3**	rectangular profile l 100 mm, 3x1 mm
	**DW-BARRA2x4**	rectangular profile l 100 mm, 4x2 mm

** 5 pieces packaging.*

Accessories for Conico system

	**PAR-PP**	parallelometer	 **PAR-INS-AS**	parallelometer insert

Remember to always purchase the parallelometer with the insert relating to the implant connection on which it is to be used.

Scanbodies for laboratory and intraoral use in titanium for Conico system	
	unique platform
	C-INT-CAMTRA-L

For further information please consult the prosthetic manual MP-IMP-PKS-ONE-E.

CAD-CAM prosthetic solutions

Customized abutments and screw retained crowns (full or with luting support)	Full bridges or bridges with luting support	DIRECT/TORONTO Full bridges or bridges with luting support	Milled bars tightened on implants
			

Scanbodies

Scanbodies for laboratory and intraoral use in titanium

	ø 3.30	ø 3.80	
	**A-INT-CAMTRA330-L**	**A-INT-CAMTRA380-L**	to transfer P.A.D. abutments connection position

Fixation screw VM2-180 included.

Scanbodies for laboratory and intraoral use in titanium for Conico system

	unique platform
	**C-INT-CAMTRA-L**

Scanbodies for laboratory and intraoral use in titanium

	unique platform	
	**PAD-INT-CAMTRA-L**	to transfer P.A.D. abutments connection position
	**A-INT-CAMTRA-PLAIN-L**	for PLAIN abutments
	**XA-INT-CAMTRA-L**	for XA posts

Scanbody for implants: fixation screw VM2-180 included. Scanbody for P.A.D. abutments: fixation screw PAD-VP-140 included. Scanbody for PLAIN and XA abutments: fixation screw A-PLAIN-VP200 included.

Supports for luting technique

Standard T-Connect

	cone height	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	4.00	**A-BASTZR-S-330-4**	**A-BASTZR-S-380-4**	**AS-BASTZR-S-425-4**	**AS-BASTZR-S-500-4**	repositionable
	6.00	**A-BASTZR-S-330-6**	**A-BASTZR-S-380-6**	**AS-BASTZR-S-425-6**	**AS-BASTZR-S-500-6**	
	4.00	**A-BASTZR-M-330-4**	**A-BASTZR-M-380-4**	**AS-BASTZR-M-425-4**	**AS-BASTZR-M-500-4**	rotating
	6.00	**A-BASTZR-M-330-6**	**A-BASTZR-M-380-6**	**AS-BASTZR-M-425-6**	**AS-BASTZR-M-500-6**	

VM2-180 fixation screw included.

T-Connect with 4.00 mm cone height and different transgingival heights

	transgingival h	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	1.00	**A-TCR-S-330-4-1**	**A-TCR-S-380-4-1**	-	-	repositionable
	2.00	**A-TCR-S-330-4-2**	**A-TCR-S-380-4-2**	-	-	
	3.00	**A-TCR-S-330-4-3**	**A-TCR-S-380-4-3**	-	-	
	1.00	**A-TCR-M-330-4-1**	**A-TCR-M-380-4-1**	-	-	rotating
	2.00	**A-TCR-M-330-4-2**	**A-TCR-M-380-4-2**	-	-	
	3.00	**A-TCR-M-330-4-3**	**A-TCR-M-380-4-3**	-	-	

VM2-180 fixation screw included.

Sleeves for wax-up modeling

	cone height	ø 3.30	ø 3.80	ø 4.25	ø 5.00
	4.00	**A-CCBAS-330-4**	**A-CCBAS-380-4**	**A-CCBAS-425-4**	**A-CCBAS-500-4**
	6.00	**A-CCBAS-330-6**	**A-CCBAS-380-6**	**A-CCBAS-425-6**	**A-CCBAS-500-6**

Interfase Dynamica

	cone height	ø 3.30	ø 3.80	
	4.00	**IND3PKR330/TIA**	**IND3PKR380/TIA**	non rotating
	4.00	**IND3PKH330/TIA**	**IND3PKH380/TIA**	repositionable

TPDH18L66 fixation screw not included.

Titanium sleeves for P.A.D. abutments

	**PAD-BASTZR**

PAD-VP-140 fixation screw included.

Supports for luting chairside technique

Chairside T-Connect supports with 5.00 mm h cementation cone					
	transgingival height	ø 3.30	ø 3.80	ø 4.25	ø 5.00
	0.50	**A-TCS-330S-C**	**A-TCS-380S-C**	**AS-TCS-425S-C**	**AS-TCS-500S-C**
	1.50	**A-TCS-330S-C-15**	**A-TCS-380S-C-15**	-	-
	2.50	**A-TCS-330S-C-25**	**A-TCS-380S-C-25**	-	-

VM2-180 fixation screw included.

Titanium premilled posts					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-CAMT-330-12-SW	**A-CAMT-380-12-SW**	**AS-CAMT-425-12-SW**	**AS-CAMT-500-12-SW**	repositionable

VM2-180 fixation screw included.

Titanium premilled for XA abutments	
	unique platform
	XA-CAMT-ABU-12-SW

A-PLAIN-VP100 fixation screw included.

Solutions for overdenture

Locator abutments and accessories

Locator abutments



transgingival h	ø 3.30	ø 3.80
1.00	**1670**	**1675**
2.00	**1671**	**1676**
3.00	**1672**	**1677**
4.00	-	**1678**
5.00	**1674**	**1679**
6.00	-	**1680**

Retainers and spacer rings for Locator abutments

	unique platform	
	**8515**	low retention capacity for impression taking, black, in polyethylene, 4 pcs packaging
	**8524**	retention 5 lb = 2268 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	**8527**	retention 3 lb = 1361 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	**8529**	retention 1,5 lb= 680 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	**8547**	retention 4 lb= 1814 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	**8915**	retention 2 lb= 907 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	**8548**	retention 1 lb= 450 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	**8514**	spacer rings in silicon rubber, 20 pcs packaging
	**8519-2**	kit: 2 Gr. 5 titanium caps, 2 spacer rings, 2 black retainers and 2 nylon retainers for each of the 3 different retention capacities
	**8540-2**	kit: 2 Gr. 5 titanium caps, 2 spacer rings, 2 black retainers and 2 nylon retainers for each of the 3 different retention capacities meant for disparallalelisms
	**8550-2**	kit: 2 Gr. 5 steel caps, 2 spacer rings, 2 black retainers and 2 nylon retainers for each of the 3 different retention capacities

Accessories for Locator abutments

	unique platform	
	**8530**	4 aluminium analogs
	**8505**	4 aluminium transfers, 4 retainers cod. 8515 included
	**8517**	4 parallelism pins in nylon
	**9530**	steel plate AISI 316 L

Locator Core Tool

	**8393**	complete Locator Core Tool
	**8397**	steel push rod for inserting the retainers in the caps
	8390	steel driver for abutment screwing/ unscrewing
	8394	retention jacket for the driver, to transfer the Locator abutments into the oral cavity

Locator R-Tx abutments and accessories

Locator R-Tx abutments



transgingival h	ø 3.30	ø 3.80
1.00	**31110-01**	**31111-01**
2.00	**31110-02**	**31111-02**
3.00	**31110-03**	**31111-03**
4.00	**31110-04**	**31111-04**
5.00	**31110-05**	**31111-05**
6.00	**31110-06**	**31111-06**

Accessories for Locator R-Tx abutment

	unique platform	
	**30013**	pack of 4 caps and 4 black retainers for impression taking
	**30001**	grey, without retention, pack of 4 pcs
	**30002**	blue, low retention, pack of 4 pcs
	**30003**	pink, medium retention, pack of 4 pcs
	**30004**	white, high retention, pack of 4 pcs
	**30012**	black, for the impression taking, pack of 4 pcs
	**08514**	spacer rings for the prosthetis relining phase, pack of 20 pcs

Accessories for Locator R-Tx abutment

	unique platform	
	**30014**	aluminium analogs for abutments ø 3.35 mm pack of 4 pcs
	**30015**	aluminium analogs for abutments ø 4.00, pack of 4 pcs
	**30016**	aluminium analogs for abutments ø 5.00, pack of 4 pcs
	**30017**	transfers with low retention, pack of 4 pcs
	30018	spacer rings for relining phase, pack of 4 pcs
	**30021**	Locator Core Tool for Locator R-Tx
	**08008**	RA screwdriver for Locator R-Tx, hex 1.25 mm
	**08007**	digital screwdriver for Locator R-Tx, hex 1.25 mm


Ball attachments and accessories

Ball attachments			
	transgingival h	ø 3.30 ●	ø 3.80 ●
	1.00	**A-AS-330-1**	**A-AS-380-1**
	2.00	**A-AS-330-2**	**A-AS-380-2**
	4.00	**A-AS-330-4**	**A-AS-380-4**

Caps for ball attachments			
		unique platform	
		**CAP-TFL-1**	polyamide cap
		**CONT-CAP-TFL-1**	steel container for polyamide cap
		**CAP-TIT-1**	Gr. 5 Ti cap with cap in two parts, Ti retention spring and plastic
		**AN-CAP-TIT-1**	spare plastic ring
		**MOL1-CAP-TIT-1**	spare steel spring, average hardness
		**MOL2-CAP-TIT-1**	spare steel spring, soft
		**CAP-1**	gold alloy cap, with 3 plastic rings and a tin spacer for the lab

Analog for ball attachments	
	unique platform
	ANAS

Bars

Overdenture on bars		
	**BARC-CAV-TIT**	castable bar, l 5.00 cm, h 3.00 mm, thickness 2.20 mm, ovoid-shaped profile with spacer
	**CAV-TIT**	divisible bar attachment in titanium for oval bars, h 3.00 mm, thickness 2.20 mm
	**BARC**	castable bar, l 5.00 cm, ø 2.20 mm
	**CAV-375**	bar attachment in gold alloy for round bars ø 2.20 mm

Clinical and laboratory prosthetic screws

Standard fixation screws		
	**VM2-180**	
	VM2-180-10	10 pcs packaging

Conical support fixation screws		
	**L-VMS-180**	
	L-VMS-180-10	10 pcs packaging

Fixation screws for Dynamic Abutment and Interfase Dynamica	
	**TPDH18L66**

Screws for Pick-up transfer		
	**VTRA2-180-15**	h 15.00
	**VTRA2-180-20**	h 20.00
	**VTRA2-180-MAN**	for manual screwing

Screws for closed-tray transfers		
	VTRA2-180-10	h 10.00

Screws for XA technique		
	**SH-VABU180-200-4**	for post h. 4.50 for screw-retained prosthesis
	SH-VABU180-200-5	for post h. 5.50 for screw-retained prosthesis
	SH-VABU180-200-6	for post h. 6.50 for screw-retained prosthesis
	A-PLAIN-VP200	for cemented prosthesis
	L-VMS-180	for cemented prosthesis

Screws for P.A.D. abutments		
	**PAD-VM-180**	
	PAD-VM-180-10	10 pcs packaging
	**PAD-VTRAL-140**	for Pick-up transfer, long
	**PAD-VTRA-140**	for Pick-up transfer, standard
	**PAD-VTRAL-140-MAN**	for Pick-up transfer, manual screwing
	**PAD-VP-140**	for prosthetic components on P.A.D. abutments
	PAD-VP-140-10	10 pcs packaging

Screws for PLAIN abutments		
	**A-PLAIN-VP200**	
	A-PLAIN-VP200-10	10 pcs packaging
	**A-PLAIN-VTRA200**	for PLAIN abutment transfers

Screws for CAD-CAM custom made prosthesis		
	**A-CAMTVABU180**	for zirconium oxide single crowns and ECHO overstructures
	**PAD-VCAM140**	for zirconium oxide ECHO overstructures on P.A.D. abutments
	**CAMPRON205-10**	spare ring washers for the head of the screw

Recommended torque for tightening prosthetic components	
8-10 Ncm	surgical cover screws, healing abutments
20-25 Ncm	prosthetic screws
25-30 Ncm	direct screwing on implants prosthetic components
8-10 Ncm	transfer screws

Warning
Fixation screws are always included in the packaging of each prosthetic component, unless otherwise indicated.
Please use laboratory screws for laboratory phases and keep the screw included with the prosthetic component for the final fastening in the patient's mouth.

For further information please consult the prosthetic manual MP-IMP-PKS-ONE-E.



FUTURE
DENTAL

FUTURE DENTAL
България, гр. София
кв. Кръстова вада
ул. „ген. Иван Колев“ 14
тел.: 088 3333 730